

ANNA LOGOR™ COSMETICS

"Classical European Care"

Any type of skin

Salon cosmetology procedure.

Classical European care is an important way to restore any type of skin. The focus of this line is cleaning and exfoliating the skin, restoring the water balance, rejuvenating and improving the complexion. To achieve the result, "Classical European Care" includes 4 important procedures: pilling, massage, restoration and rejuvenation. As a consequence, after the procedure, the skin is restored, looks healthy and rejuvenated.

Sequence of the procedure:

1. Prepare the skin by thoroughly cleansing with **Foam Facial Cleanser**. Remove makeup from the eyes with a damp cotton pad using **Make-up Remover**.
2. Wipe your face with **Purifying Complex Toner**.
3. In a circular motion, apply **Massage Gel** (Oily, Normal /Oily skin) or **Massage Cream** and massage for 5-10 minutes.
4. To exfoliate the skin, apply **Micro-dermabrasion Cream** and leave for 1 - 2 minutes. Starting from the forehead and moving down, holding the head with one hand, gently rub it with the other hand, creating a slight tension. Remove as many products as possible. Remove the remainder with a damp sponge.
5. Apply **Pre-steam Gel** 1-2 mm thick. Before applying, protect eyes with a cotton pad and leave on the face for 7-10 minutes. Conduct a manual massage (optional).
6. Carry out a thorough cleaning of your face with a damp towel or cotton pad.
7. Apply a moderate amount (1-2 mm thick) of the **Purifying Complex Mask** for oily and normal skin or **Cranberry Mask** for dry and fading skin and leave for 10 minutes.
8. Remove the mask with a warm, damp cloth.
9. Wipe your face with **Purifying Complex Toner**.
10. Apply **Vitamin C Serum**. To improve efficiency, use the High Frequency Current device for 2 -3 minutes.
11. Apply a small amount of **Copper Plus Moisturizer** (face and neck area) and **Daily Eye Therapy Cream**(eye area). With the help of circular, massage movements, fully incorporate the cream into the skin of the face and neck.
12. To get the procedure, apply a sufficient amount of **Mezotherapy Gel** on the face and neck. Conduct manual massage for up to 3 minutes to facilitate the penetration of active ingredients into the skin.

Cosmetic products necessary for the procedure

"Classical European care":

1. Purifying Complex Toner	(ART# 302)	7. Regeneration Eye Cream	(ART# 722)
2. Purifying Complex Mask	(ART# 355)	8. Micro-Dermabrasion Cream	(ART# 924)
3. Copper Plus Moisturizer	(ART# 423)	9. Make-up Remover	(ART# 1001)
4. Silky Cranberry Mask	(ART# 455)	10. Pre-Steam Gel	(ART# 1002)
5. Mezotherapy Gel	(ART# 582)	11. Foam Facial Cleanser	(ART# 1003)
6. Vitamin C Serum	(ART# 653)	12. Massage Cream	(ART# 1083)
		13. Massage Gel	(ART# 1084)