

ANNA LOGOR™ COSMETICS

«Ultra-moisturizing»

Any type of skin
Salon cosmetology procedure.

Water is an essential element for the skin, so that it looks elastic and healthy. Maintaining a natural water balance is a must. Dehydrated skin loses its elasticity, which contributes to the appearance of wrinkles. This moisturizing treatment for the face is an intensive method of restoring the water balance, as a result of which the skin of the face is moisturized and as a result - the tone and elasticity improve, the skin becomes tender, velvety and acquires a healthy appearance.

Sequence of the procedure:

1. Prepare the skin by thoroughly cleansing the face of **Foam Facial Cleanser**. Remove the make-up around your eyes with a damp cotton pad with **Make-up Remover**.
2. Apply **Shuttle Spray**, neatly sprayed on the face, in order to prepare the skin for further procedures.
3. Perform a manual massage for relaxation for 5-10 minutes, using **Massage Gel** (oily, combination skin) or **Massage Cream**.
4. Protect your eyes with cotton pads and apply a layer 1-2 mm thick **Pre-Steam Gel**. Steam the face for 7-10 minutes. Manual massage - if desired.
5. Thoroughly clean the face of the applied gel with a warm, damp towel.
6. Re-apply the **Shuttle Spray**.
7. Apply **Hydrating Serum** and incorporate into the skin using an ULTRASONIC device within 5 to 7 minutes. (Manual massage for 2-3 minutes can be used in the absence of an Ultrasound machine).
8. Apply a moderate amount of **Hydrating Mask** and leave for 10 minutes.
9. Remove the mask with a warm, damp cloth or cotton pad.
10. Re-apply the **Shuttle Spray**.
11. Using a High Frequency Current device, lightly massage for 2-3 minutes.
12. Apply a small amount of **Hydrating Moisturizer** (on the face skin) and **Regeneration eye cream** (eye area). With the help of circular, massage movements fully incorporate the cream into the skin of the face and neck.
13. To complete the procedure, apply a sufficient amount of **Mezotherapy Gel** on the face and neck. Conduct manual massage for up to 3 minutes to facilitate the penetration of active ingredients into the skin.

Cosmetic products necessary for the procedure

" Ultra-moisturizing":

1. Mezotherapy Gel	(ART# 582)	7. Make-up Remover	(ART# 1001)
2. Shuttle Spray	(ART# 602)	8. Pre-Steam Gel	(ART# 1002)
3. Hydrating Moisturizer	(ART# 621)	9. Foam Facial Cleanser	(ART# 1003)
4. Hydrating Serum	(ART# 633)	10. Massage Cream	(ART# 1083)
5. Hydrating Mask	(ART# 645)	11. Massage Gel	(ART# 1084)
6. Regeneration eye cream	(ART# 722)		